OSTENIL®

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Batch number

OSTENIL®

Sodium hyaluronate 1.0%. Viscoelastic solution for injection into the joint cavity. Sterile by moist heat.

Composition:

Each pre-filled syringe of 20 mg/2.0 mL contains:

Principal component: sodium hyaluronate 1.0%

Other components: sodium chloride, disodium phosphate, sodium dihydrogen phosphate, water for injections.

Indications:

Temporary synovial fluid replacement for lubricating in degenerative and traumatic changes of knee joint and other synovial joints resulting in pain decrease and joint mobility improvement.

Contra-indications:

OSTENIL® should not be used in patients with ascertained hypersensitivity to any of its constituents.

Interactions:

No information on the incompatibility of OSTENIL® with other solutions for intra-articular use is available to date. The concomitant use of an oral analgesic or anti-inflammatory drug during the first few days of treatment may be helpful for the patient.

Undesirable effects:

Local secondary phenomena such as pain, feeling of heat, redness and swelling may occur in the joint treated with OSTENIL®. Application of an ice pack for five to ten minutes onto the treated joint will reduce the incidence of these events.

Dosage and administration:

The clinically recommended dose is 3 injections of 2.0 mL (20 mg) of OSTENIL® into the affected joint. Several joints can be treated at the same time. Depend on situation and symptoms of disease, doctor can request up to 5 consecutive doses at 1 week interval. Repeated treatment cycles may be administered as required.

In case of joint effusion it is advisable to reduce the effusion by aspiration, rest, application of an ice pack and/or intra-articular corticosteroid injection. Treatment with OSTENIL® can be started two to three days later.

The content and the outer surface of the OSTENIL® pre-filled syringe are sterile as long as the sterile pack is intact. Take the pre-filled syringe out of the sterile pack, unscrew the Luer lock cap from the syringe, attach a suitable needle (for example 19 to 21 G) and secure it by turning slightly. Remove any air bubble, if present, before injection.

Precautions:

Caution should be exercised in patients with known hypersensitivity to one of the constituents. As with all invasive joint treatments in very rare cases an infection may occur. Hence, the general precautions for intra-articular injections should be observed, including measures to avoid joint infections. OSTENIL® should be injected accurately into the joint cavity, if necessary under imaging control. Avoid injections into blood vessels or surrounding tissues! As no clinical evidence is available on the use of hyaluronic acid in children, pregnant and lactating women or in inflammatory joint diseases such as rheumatoid arthritis or Bechterew disease, treatment with OSTENIL® is not recommended in these cases. Do not use if the pre-filled syringe or sterile pack are damaged. Any solution not used immediately after opening must be discarded. Otherwise the sterility is no longer guaranteed and this may be associated with a risk of infection. Store below 30 °C! Do not use after the expiry date indicated on the box. Keep out of the reach of children.

Characteristics and mode of action:

Synovial fluid, which is viscoelastic due to the presence of hyaluronic acid, is found in all synovial joints, particularly the large weight bearing joints, where it ensures normal, painless movement due to its lubricating and shock-absorbing properties. It is also responsible for the nutrition of the cartilage. In degenerative joint disorders such as osteoarthritis, the viscoelasticity of the synovial fluid is markedly reduced thereby decreasing its lubricating and shock-absorbing functions. This increases mechanical loading of the joint and cartilage destruction which ultimately results in pain and restricted mobility of the affected joint. Supplementing this synovial fluid with intra-articular injections of highly purified hyaluronic acid can ameliorate the viscoelastic properties of synovial fluid. This improves its lubricating and shock-absorbing functions and reduces mechanical overload of the joint. As a rule this results in a decrease in pain and an improvement in joint mobility which may last for several months after a treatment cycle of five intra-articular injections.

Presentation:

One pre-filled syringe of 20 mg/2.0 mL sodium hyaluronate in a sterile pack.

Three pre-filled syringes of 20 mg/2.0 mL sodium hyaluronate in sterile packs.

Five pre-filled syringes of 20 mg/2.0 mL sodium hyaluronate in sterile packs.

OSTENIL® is a medical device. To be used by a physician only.

Last revision date: 2022-07

Refer to instructions for use

For single use only

Expiry date

Sterile by moist heat

Don't use if the sterile barrier
is damaged

Store below 30 °C!

Product owner:

TRB CHEMEDICA AG
26 Otto-Lilienthal-Ring
85622 Feldkirchen (München), Germany

Manufacturer:

TRB CHEMEDICA SA
Chemin St-Marc 3,
1696 Vouiry, Switzerland

Origin: Switzerland

Importer and Distributor:

DKSH Vietnam Co., Ltd.